






3R78

ZH 使用说明书	4
DE Gebrauchsanweisung	10
EN Instructions for Use	18
FR Instructions d'utilisation	25
IT Istruzioni per l'uso	33
ES Instrucciones de uso	42
PT Manual de utilização	50
NL Gebruiksaanwijzing	58
SV Bruksanvisning	66
DA Brugsanvisning	73
NO Bruksanvisning	80
FI Käyttöohje	88
PL Instrukcja użytkowania	95
HU Használati utasítás	103
CS Návod k použití	111
RO Instrucțiuni de utilizare	118
HR Upute za uporabu	126
TR Kullanma talimatı	134
JA 取扱説明書	141
KO 사용 설명서	148

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CAUTION	Warnings regarding possible risks of accident or injury.
NOTICE	Warnings regarding possible technical damage.
INFORMATION	Additional information on the fitting/ use.

INFORMATION Date of the last update: 2014-04-16 <ul style="list-style-type: none"> • Please read this document carefully. • Follow the safety instructions.
--


1 Description

1.1 Intended Use

The 3R78 polycentric knee joint is intended exclusively for use in prosthetic fittings of the lower extremities.

1.2 Area of Application

Area of application according to the Ottobock MOBIS mobility system:

Recommended for amputees with **mobility grades 2 and 3**
(restricted outdoor walkers and unrestricted outdoor walkers).

Approved for a patient weight of up to 100 kg

1.3 Environmental Conditions

Allowable environmental conditions
Temperature range for use: -10 °C to +45 °C
Allowable relative humidity 0 % to 96 %, non-condensing

Unallowable environmental conditions
Mechanical vibrations or impacts
Perspiration, urine, fresh water, salt water, acids
Dust, sand, highly hygroscopic particles (e. g. talcum)

1.4 Safety Instructions

Please make your patients aware of the following safety instructions:

CAUTION Risk of falling as a result of excessive strain. The Ottobock prosthetic components described here have been developed for prosthesis wearers with mobility grades 2 and 3 and must not be used by prosthesis wearers with higher mobility grades. Otherwise, malfunctions of the joint can occur which may result in breakage of load-bearing parts. This can cause the patient to fall. Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures the patient's safety.
--

Should prosthetic components be subjected to extreme strain (such as a fall), contact your prosthetist immediately and have the components inspected for possible damage. The responsible prosthetist is your point of contact and will pass the prosthesis on to Ottobock Service if applicable.

CAUTION

Risk of falling due to incorrect environmental conditions. Using the knee joint in incorrect environmental conditions can lead to malfunctions of the joint up to loss of functionality due to structural failure. This can cause the patient to fall.

Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids. Penetration of liquids, dust, sand, etc. as well as the use of the knee joint at places with excessive ambient moisture must be avoided.

NOTICE

Risk of damage due to incorrect environmental conditions. Prosthetic components must not be exposed to corrosive elements such as fresh water, salt water and acids. Penetration of liquids, dust, sand, etc. as well as the use of the knee joint at places with excessive ambient moisture must be avoided.

Using this medical device in such environmental conditions will render all claims against Otto Bock HealthCare null and void.

1.5 Construction and Function

The 3R78 prosthetic knee joint is a polycentric knee joint with pneumatic swing phase control and integrated extension assist spring.

Joint structure

The proximal joint connection is realised by means of a pyramid adapter (Figure 1, Item A). The upper joint section (Figure 1, Item B) and lower joint section (Figure 1, Item C) are connected by the anterior links (Figure 1, Item D) and the posterior link (Figure 1, Item E) to form a four-bar linkage. In extension, the instantaneous centre of rotation is located above the joint and behind the load line, thus providing knee stability in the stance phase (Figure 1, Item F).

The distal portion of the knee joint is fabricated as a clamp fitting to accept the tube adapter (Figure 1, Item G).

Swing phase control

To control swing phase, the pneumatic unit is linked to the upper joint section with a connecting rod (Figure 2, Item A). An integrated extension assist spring provides support in this phase.

The flexion (Figure 2, Item B) and extension resistance (Figure 2, Item C) can be separately adjusted in order to meet the individual needs of the prosthesis wearer (see Section 3.2.2). Movement resistance dependent on the walking speed prevents excessive swing-through of the lower prosthetic leg in flexion and assures harmonious, dampened extension.

1.6 Combination Possibilities

Depending on the patient's functional demands, the following prosthetic feet and tube adapters can be used:

- Dynamic Motion (1D35), Greissinger plus (1A30), Trias (1C30), C-Walk (1C40), Axion (1E56); Axion DP (1E58), Adjust (1M10).
- 2R49 tube adapter (aluminium), 2R3 tube adapter (stainless steel), 2R38 tube adapter (titanium).

2 Technical Data

Order no.	3R78
Proximal connection	Pyramid adapter
Distal connection	Tube clamp, 30 mm Ø
Knee flexion angle (Figure 3)	150°
Weight	760 g
System height	156 mm
Proximal system height up to alignment reference point	-7 mm
Distal system height up to alignment reference point	163 mm
Max. user weight	100 kg / 220 lbs
Mobility grade	2, 3

3 Handling

3.1 Alignment

CAUTION

Risk of falling as a result of incorrect alignment. Errors during the alignment of the prosthesis can lead to malfunction of the joint up to loss of functionality due to structural failure. This can cause the patient to fall. Observe the alignment instructions.

The three-dimensional arrangement of the prosthetic socket and the modular components (knee joint, prosthetic foot etc.) affects the static and dynamic functions of the prosthesis.

The advantages of the 3R78 can only be optimised with correct alignment.

The optimal residual limb position must be anticipated when positioning the socket connector (Figure 6). Plumb lines in the frontal and sagittal planes, drawn from the hip joint centre of rotation and marked while taking the plaster cast and during trial fitting of the test socket, will facilitate correct positioning of the lamination anchor or socket adapter.

To align the prosthesis please proceed in 2 steps:

Step one - bench alignment on the alignment apparatus (e.g. 743L200 LASAR Assembly, 743A200 PROS.A Assembly).

Step two - alignment optimisation using the 743L100 LASAR Posture.

3.1.1 Bench Alignment with Alignment Apparatus (Figure 6)

- 1 Position the middle of the foot (MF) **approx. 30 mm anterior** to the alignment reference line. This applies to all foot components that are recommended for use with 3R78, independently of the specifications in the Instruction for Use for said feet!
- 2 Add **5 mm** to the effective heel height. Set correct outward rotation of the foot.
- 3 Clamp the knee joint using the appropriate adapter inserts. For bench alignment, the alignment reference line should run through the upper anterior axis (alignment reference point). Take into account the knee-ground distance and outward rotation of the knee (the adapter insert provides for a rotation of approx. 5°). Recommended sagittal positioning of the alignment reference point: **20 mm** above the medial tibial plateau.
- 4 Connect the prosthetic foot and modular knee joint. The connection to the prosthetic foot is realised using the recommended tube adapter.

CAUTION

Risk of falling due to incorrectly mounted tube adapter. Incorrect mounting of the tube adapter in the tube clamp can lead to excessive load on the tube clamp and to insufficient force transmission and as a result cause the patient to fall.

An insufficiently secured connection between the tube adapter and prosthetic foot will also lead to a risk of falling.

Use the recommended tube adapter (see Section 1.5).

Observe the following instructions.

- a. Do not clamp the tube into a vice! Cut the tube using 719R2 Tube Cutter placed at a right angle to the tube.

NOTICE

Risk of damage to the prosthesis: For reasons of structural strength, the residual pieces of cut tube adapters must not be used as components in prostheses.

- b. Afterwards carefully deburr the inside and outside of the tube using the 718R1 Tube Deburring Knife.
- c. Do not make any length adjustment by varying the insertion depth. Insert the tube into the knee joints' tube clamp all the way to the stop.
- d. Tighten the cap screw of the tube clamp using the 710D1 Torque Wrench. Torque: **10 Nm** (Fig. 4).
- e. Connect the tube adapter to the foot. To do so, tighten the 506G3 Set Screws which are located next to each other at the distal end of the tube adapter using the 710D4 Torque Wrench. Torque: **10 Nm**. (For completion of the prosthesis, apply 636K13 Loctite to the set screws and tighten all set screws with **15 Nm**).
- 5 Mark the centre of the socket proximally and distally on the lateral side. Draw a line through both marks from socket brim to the distal end of the socket.
- 6 Now position the socket such that the alignment reference line passes through the proximal centre mark. Set the socket flexion to somewhere between 3° and 5° taking the individual situation (e.g. hip flexion contractures) and the ischial tuberosity to ground distance into account.

CAUTION

Risk of falling as a result of incorrect alignment. If residual limb flexion is not taken into account, the joint will be positioned too far to the front. This can lead to malfunction and premature wear and cause the patient to fall.

- 7 Connect the socket and modular knee joint using a corresponding adapter (e.g. 4R111, 4R41, 4R55, 4R51 Socket Adapter).

3.1.2 Static Alignment Optimization with 743L100 L.A.S.A.R. Posture (Fig. 7)

After bench alignment of the prosthesis, static alignment using the L.A.S.A.R. Posture is carried out. In order to ensure appropriate stability combined with easy initiation of the swing phase, please proceed as follows:

- 1 To make the load line visible, the transfemoral amputee (with shoes) stands on the L.A.S.A.R. Posture with the prosthetic side on the force plate and with the other leg on the height compensation panel. The prosthetic side should be sufficiently loaded (> 35% of the body weight).

- 2 Now adapt the alignment by only adjusting the plantar flexion so that the load line (laser line) is **35 mm** anterior to the alignment reference point (= knee axis) (see Fig. 7).

3.1.3 Dynamic Alignment Optimization

Perform dynamic optimisation during trial walking. Often, the following aspects have to be observed and adapted, if necessary:

- Socket flexion position by verifying step length symmetry (sagittal plane)
- Adduction position of the socket and M-L positioning of the socket adapter (frontal plane)
- Rotation position of the knee axis and outward rotation of the prosthetic foot (transversal plane)

Finally, document the fitting result using for example the “Dartfish” motion analysis software.

3.2 Adjustment

CAUTION

Risk of pinching. When adjusting/assembling the knee joint, there is a risk of pinching in the flexion range area of the knee joint. **Do not reach into the flexion range area.**

CAUTION

Risk of falling during trial walking. For safety reasons, the first steps of the prosthesis wearer should always be taken between parallel bars.

Trial walking should initially begin with the factory default settings and the functions should be explained to the patient.

3.2.1 Stance Phase Stability

CAUTION

Risk of falling due to insufficient posterior positioning of the knee joint. The position of the knee joint in the prosthesis has to ensure that the polycentric structure always assures sufficient stance phase stability during heel impact and subsequent load transfer.

Contrary to single-axis knee joints, polycentric knee joints become stable upon heel impact strike because of their axis geometry.

Knee stability is affected by the alignment, foot construction, residual limb characteristics, joint position and location of the instantaneous centre of rotation.

When required, stance phase stability can be corrected as follows:

Situation	Measure
Movement can only be initiated with a great deal of effort	Check the alignment of the prosthesis. Check whether existing residual limb contracture is being taken into account and adapt socket flexion if applicable. Reduce plantar flexion if applicable.
Patient experiences buckling	Check the alignment of the prosthesis. Check whether existing residual limb contracture is being taken into account and adapt socket flexion if applicable. Gradually increase plantar flexion if applicable.

3.2.2 Flexion and Extension Adjustment

CAUTION

Risk of falling due to incorrectly adjusted extension resistance. When adjusting extension resistance, ensure that full extension is consistently achieved during heel impact. This is the only way to ensure the joint offers sufficient stance phase stability during subsequent load transfer.

The resistances for flexion (F) and extension (E) are individually and independently adjustable. Movement resistance dependent on the walking speed prevents excessive swing-through of the lower prosthetic leg in flexion and assures harmonious, dampened extension.

Upon delivery, the valve screw (Figure 5, Item B) for extension damping is turned out as far as possible (set to the lowest resistance) and the valve screw (Figure 5, Item A) for flexion damping is closed by 2 turns. Use a 2 mm Allen wrench (not included in the scope of delivery!) to change the setting of the valve screws.

Proceed as follows:

1. During trial walking, first adjust flexion resistance (F) (Figure 5, Item A). Approach the optimum for various walking speeds in small adjustment increments.
2. Then adjust extension resistance (E) (Figure 5, Item B). Approach the optimum in small adjustment increments.
3. Finally, fine-tune the coordination between flexion and extension.

Problem	Cause	Measure
F > flexion valve (see Figure 5, Item A)		
Lower prosthetic leg swings back too far	Flexion resistance set too low	Increase flexion resistance Turn valve screw clockwise
Lower prosthetic leg does not swing back far enough	Flexion resistance set too high	Reduce flexion resistance Turn valve screw counter-clockwise
E > extension valve (see Figure 5, Item B)		
Extension impact too hard	Extension resistance set too low	Increase extension resistance Turn valve screw clockwise
Extension movement too slow; extension stop (full extension) not reached upon heel impact	Extension resistance set too high	Reduce extension resistance Turn valve screw counter-clockwise

3.3 Final Assembly

The 3S107 Foam Cover can be used for the 3R78.

The foam cover may affect knee joint functionality. Readjustment might therefore become necessary to compensate for such effects.

CAUTION

Risk of falling when using talcum powder. Talcum powder absorbs grease from mechanical components. This leads to mechanical faults - e.g. a stiff pneumatic unit - and may cause the patient to fall.

Do not use talcum powder. Using this medical device after application of talcum powder will render all claims null and void.

INFORMATION

An alternative for eliminating noises is a silicone spray (519L5) that can be sprayed directly on the friction surface of the foam cover.

3.4 Maintenance Instructions

Please make your patients aware of the following safety instructions:

⚠ CAUTION

Risk of injury as a result of excessive strain. The knee joint function may be affected by both environmental and service conditions. To reduce the risk for the patient, the knee joint should not be used when noticeable functional changes occur.

Noticeable functional changes may include poor response, incomplete extension, diminished stance phase stability, deteriorating swing phase behaviour, abnormal noises, etc.

Measure: Visit an authorised specialist workshop to have the prosthesis inspected.

NOTICE

Risk of damage caused by improper cleaning agents. The use of aggressive cleaning agents must be avoided. They may cause damage to bearings, seals and plastic parts.

Do not disassemble the knee joint! If you have a problem, please send the joint to Ottobock.

INFORMATION

This component has been tested according to the ISO 10328 standard for three million load cycles. Depending on the amputee's activity this corresponds to a service life of three to five years.

We recommend carrying out regular safety checks once a year.

INFORMATION

When using exoprosthetic knee joints, servomotoric, hydraulic, pneumatic or brake load dependent control functions can cause movement noise. This kind of noise is normal and unavoidable. Usually, it does not cause any problems.

If movement noise increases noticeably during the life cycle of the the knee joint, the joint should be inspected by a prosthetist immediately.

Ottobock recommends readjusting the knee joint settings once the patient has gotten used to the prosthesis. This period of time varies depending on individual patient characteristics.

Please check the knee joint for wear and proper function at least once a year and make adjustments if necessary. Special attention should be paid to movement resistance and abnormal noises. Full flexion and extension must be guaranteed.

Do not disassemble the joint. If there are any potential malfunctions, please send the entire joint to Ottobock Service.

4 Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The device must only be operated with tested modular components in accordance with the Ottobock MOBIS mobility system. The manufacturer is not liable for damage caused by component combinations that were not authorised by the manufacturer.

5 CE Conformity

This device meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Ottobock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH